



単位 mm

| A | B | C | D | E | F | G |
|------------|-----------|-----------|-----------|-----------|-----------|-----------|
| 110.0 ±0.3 | 30.0 ±0.2 | 20.0 ±0.2 | 18.0 ±0.2 | 16.0 ±0.2 | 14.5 ±0.2 | 25.5 ±0.2 |

| | | | |
|-------------|------|------|----------------|
| コードNo. | 4088 | 品名 | PP分配管 8mm x 5連 |
| 2004年11月12日 | 版 1 | 株式会社 | サンプラテック |