




(mm)

| A    |      | B    |      | C    |      | D    |      | E    |      | F   |      | G    |      |
|------|------|------|------|------|------|------|------|------|------|-----|------|------|------|
| 92.0 | ±1.0 | 47.9 | ±0.5 | 32.8 | ±0.3 | 21.8 | ±0.2 | 46.2 | ±0.2 | 8.0 | ±0.2 | 36.0 | ±0.2 |

code No.

18112

サンプラ<sup>®</sup> PFA 広口ボトル (中栓付) 100ml  
PFA WIDE MOUTH BOTTLE(with inner seal) 100ml

13th November,2024

Ver.2

株式会社 **サンプラテック**  
 **SANPLATEC CORP.**